

SPRITE

DIRECT CONNECT HPLC COLUMN



Guide to SPRITE Part Numbers

xR-0121-xxxx 10 x 2.1mm Sprite
 xR-0221-xxxx 20 x 2.1mm Sprite
 xR-0421-xxxx 40 x 2.1mm Sprite

See Page 23 in the Higgins Catalog for complete Part Number information

The Sprite family of narrowbore HPLC columns is a unique high performance product for LC-MS and other high-speed, small volume applications. All three sizes of 2.1mm I.D Sprite HPLC columns are individually tested. A Certificate of Analysis accompanies each column. The Sprite's certificate includes the original chromatogram and comprehensive testing and sorbent documentation.

The inert PEEK and titanium design has been proven to be structurally sound to over 5,000 PSI (340Bar).

The optional Higgins *Micro Column Connection Kit* (P/N MIC-KIT, see page 69 in the Higgins Analytical Catalog) includes sleeves and finger-tight bushings to facilitate the attachment of a Sprite to most injection valves, detectors, and MS source inlets.

CLIQUEUS Sprite



CLIQUEUS C8 5µm

10X2.1mm CR-0121-C085	\$190	20X2.1mm CR-0221-C085	\$200	40X2.1mm CR-0421-C085	\$235
--------------------------	-------	--------------------------	-------	--------------------------	-------

CLIQUEUS C18 5µm

10X2.1mm CR-0121-C185	\$190	20X2.1mm CR-0221-C185	\$200	40X2.1mm CR-0421-C185	\$235
--------------------------	-------	--------------------------	-------	--------------------------	-------

CLIQUEUS Cyano 5µm

10X2.1mm CR-0121-CNP5	\$190	20X2.1mm CR-0221-CNP5	\$200	40X2.1mm CR-0421-CNP5	\$235
--------------------------	-------	--------------------------	-------	--------------------------	-------

CLIQUEUS Phenyl 5µm

10X2.1mm CR-0121-PHN5	\$190	20X2.1mm CR-0221-PHN5	\$200	40X2.1mm CR-0421-PHN5	\$235
--------------------------	-------	--------------------------	-------	--------------------------	-------

CLIQUEUS Silica 5µm

10X2.1mm CR-0121-SIL5	\$190	20X2.1mm CR-0221-SIL5	\$200	40X2.1mm CR-0421-SIL5	\$235
--------------------------	-------	--------------------------	-------	--------------------------	-------

HAISIL 100 Sprite



HAISIL 100 C8 5µm

10X2.1mm HR-0121-M085	\$190	20X2.1mm HR-0221-M085	\$200	40X2.1mm HR-0421-M085	\$235
--------------------------	-------	--------------------------	-------	--------------------------	-------

HAISIL 100 C18 5µm

10X2.1mm HR-0121-M185	\$190	20X2.1mm HR-0221-M185	\$200	40X2.1mm HR-0421-M185	\$235
--------------------------	-------	--------------------------	-------	--------------------------	-------

HAISIL 300 Sprite



HAISIL 300 C18 5µm

10X2.1mm HR-0121-W185	\$190	20X2.1mm HR-0221-W185	\$200	40X2.1mm HR-0421-W185	\$235
--------------------------	-------	--------------------------	-------	--------------------------	-------



HAISIL HL Sprite



HAISIL HL C18 5µm

10X2.1mm LR-0121-C185	\$190	20X2.1mm LR-0221-C185	\$200	40X2.1mm LR-0421-C185	\$235
--------------------------	-------	--------------------------	-------	--------------------------	-------

HEAVY Sprite



HEAVY C18 5µm

10X2.1mm VR-0121-C185	\$190	20X2.1mm VR-0221-C185	\$200	40X2.1mm VR-0421-C185	\$235
--------------------------	-------	--------------------------	-------	--------------------------	-------

PHALANX Sprite



PHALANX C18 5µm

10X2.1mm PR-0121-C185	\$190	20X2.1mm PR-0221-C185	\$200	40X2.1mm PR-0421-C185	\$235
--------------------------	-------	--------------------------	-------	--------------------------	-------

PROTO 200 Sprite

200

PROTO 200 C4 5µm

10X2.1mm RR-0121-D045	\$190	20X2.1mm RR-0221-D045	\$200	40X2.1mm RR-0421-D045	\$235
--------------------------	-------	--------------------------	-------	--------------------------	-------

PROTO 200 C18 5µm

10X2.1mm RR-0121-D185	\$190	20X2.1mm RR-0221-D185	\$200	40X2.1mm RR-0421-D185	\$235
--------------------------	-------	--------------------------	-------	--------------------------	-------

PROTO 300 Sprite



PROTO 300 C4 5µm

10X2.1mm RR-0121-W045	\$190	20X2.1mm RR-0221-W045	\$200	40X2.1mm RR-0421-W045	\$235
--------------------------	-------	--------------------------	-------	--------------------------	-------

PROTO 300 C18 5µm

10X2.1mm RR-0121-W185	\$190	20X2.1mm RR-0221-W185	\$200	40X2.1mm RR-0421-W185	\$235
--------------------------	-------	--------------------------	-------	--------------------------	-------

TARGA Sprite



TARGA C8 5µm

10X2.1mm TR-0121-C085	\$190	20X2.1mm TR-0221-C085	\$200	40X2.1mm TR-0421-C085	\$235
--------------------------	-------	--------------------------	-------	--------------------------	-------

TARGA C18 5µm

10X2.1mm TR-0121-C185	\$190	20X2.1mm TR-0221-C185	\$200	40X2.1mm TR-0421-C185	\$235
--------------------------	-------	--------------------------	-------	--------------------------	-------



SPRITE is Designed to Operate in Either Flow Direction

A SPRITE is made of PEEK and titanium and can be connected directly to a sample injection valve or directly to a mass spectrometer's source inlet.

