



# Capillary™

*A New Benchmark for Capillary Column  
Quality and Cost*

- ✓ All 75µm and 150µm columns are 250mm long
- ✓ Use as is or customize with included ceramic cutter
- ✓ Save operating expenses - costs a fraction of competitors'
- ✓ Reduce variability with phases renowned for their reliable stability and lot-to-lot reproducibility

- ✓ Enormous selection of chemistries

Twelve polymeric RP and ion exchangers

Over 20 silica based sorbents including:

**TARGA, CLIPEUS, PROTO200, PROTO300, PHALANX**

TARGA C18 5µm  
250 x 0.075mm  
P/N: TC-2507-C185  
S/N 166435



Higgins  
Analytical