



Capillary™

*A New Benchmark for Capillary Column
Quality and Cost*

- All 75µm and 150µm columns are 250mm long
- Use as is or customize with included ceramic cutter
- Save operating expenses - costs a fraction of competitors'
- Reduce variability with phases renowned for their reliable stability and lot-to-lot reproducibility
- Enormous selection of chemistries

*Twelve polymeric reverse phases, ion exchangers and
Over 20 silica based sorbents including:*

TARGA, CLIPEUS, PROTO200, PROTO300, PHALANX

TARGA C18 5µm
250 x 0.075mm
P/N: TC-2507-C185
S/N 166435



Higgins
Analytical